



PROMETHEUS  
MEDICAL LTD (MENA)

# Blast™ Bandage

The Blast™ Bandage provides the ability to quickly package traumatic amputations, burns and large pattern wounds with minimal effort. The Blast™ Bandage provides a 47cm x 47cm treatment area, but packs to about the size of a 10cm combat bandage.

Following traumatic amputation, the large non-adherent wound pad is used to gather together the injured tissue of the stump, before pressure is applied using the integral elastic bandage. Covering and applying effective haemostatic stump -pressure is often very challenging. The Blast™ Bandage is specifically designed for this problem, making the task quick and simple to achieve.



The Blast™ Bandage also features a large removable occlusive dressing. This can be used to cover eviscerated bowel (keeping it moist, warm and clean); as a burns dressing (being more robust than cling film); and an occlusive dressing for multi-penetrating chest wounds.

The 'brakes' on the Blast™ Bandage help prevent accidental unrolling of the elastic wrap during application. These brakes are especially helpful during the dressing of amputations. By providing multiple points of adhesion the brakes prevent slippage of the elastic during the wrapping process. This makes the bandaging process much easier to achieve, with a more secure end-result.

Specifications	
Size Packed	10cm x 10cm
Size Unpacked	47cm x 47cm
Weight	100g
Price	Available on request

Prometheus Medical Ltd  
The Old Rectory, Hope Under Dinmore,  
Herefordshire HR6 0PW  
T: +44(0)1568 613942 F: +44(0)1568 620032

Prometheus Medical Ltd (MENA)  
1st Floor, Al Ahli Bank Building, Al Wakra Branch,  
Al Wakra Main Street, PO Box 172222, Doha, Qatar  
T: +974 447 700 73/5 F: +974 447 700 71